

## Product Datasheet

### Transferrin Antibody JP148

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Transferrin Antibody                                       |
| <b>Catalog Number</b>      | JP148  |
| <b>Species/Host</b>        | Goat   |
| <b>Tested Applications</b> | TIA  |
| <b>Immunogen</b>           | Purified Human Transferrin                                 |
| <b>Brand</b>               | EastCoast Bio  |
| <b>Form/Appearance</b>     | Polyclonal, Delipidated and 0.2 µm filtered                |
| <b>Storage</b>             | 2-8°C  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Buffer</b>              | 20 mM Sodium Phosphate, pH 7.2 with 150 mM Sodium Chloride |
| <b>Specificity</b>         | Single precipitin arc visible with immunoelectrophoresis   |
| <b>Preservatives</b>       | 0.05% Sodium Azide   |