

Product Datasheet

Poliovirus-HRP Antibody JP149

Product Name	Poliovirus-HRP Antibody
Catalog Number	JP149
Species/Host	Goat
Tested Applications	ELISA, WB, IHC
Immunogen	Purified Poliovirus types 1-3
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody-HRP Conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	10 mM Phosphate Bufered Saline, pH 7.2 with 10 mg/ml BSA
Specificity	Poliovirus Types 1-3
Preservatives	0.002% Thimerosal