

Product Datasheet

Apolipoprotein CIII (Apo CIII) Antibody JP184

Product Name	Apolipoprotein CIII (Apo CIII) Antibody
Catalog Number	JP184
Description	Affinity Purified Polyclonal Antibody
Species/Host	Goat
Tested Applications	WB, ELISA
Immunogen	Human Apo CIII
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	75 mM Sodium Phosphate, 75 mM NaCl, 0.5 mM EDTA, pH 7.2
Specificity	Human Apo CIII
Preservatives	0.02% Sodium Azide