

Product Datasheet

Borrelia burgdorferi Antibody, FITC JP191

Product Name	Borrelia burgdorferi Antibody, FITC
Catalog Number	JP191
Species/Host	Rabbit
Tested Applications	IFA
Immunogen	B. burgdorferi whole cells
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, Fluorescein Conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 with BSA
Cross Reactivity	Cross reacts with Treponema pallidum, B. hermsii, and B. parkerii
Specificity	Reactive with WB bands present at 83kD, 41kD, 34kD, 31kD and other low molecular weight bands.
Preservatives	0.1% Sodium Azide