

# Product Datasheet

## Human IL-1 beta ELISA Set 851.610.010

Product Name	Human IL-1 beta ELISA Set
Catalog Number	851.610.010
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	15.6-500 pg/ml
Sensitivity	7 pg/ml
Tested Applications	ELISA
Alternative Names	IL-1b, IL-1beta
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural (pro and mature) and recombinant human IL-1b
Sample Types	Serum, Plasma, Cell culture supernatant