

Product Datasheet

Transferrin Antibody JP202

Product Name	Transferrin Antibody
Catalog Number	JP202
Species/Host	Goat
Tested Applications	TIA
Immunogen	Purified Human Transferrin
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	50 mM Tris, pH 7.4 with 500 mM NaCl
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.09% Sodium Azide