
Product Datasheet

Chloramphenicol-BSA Antigen LA005

Product Name	Chloramphenicol-BSA Antigen
Catalog Number	LA005
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	0.1M Phosphate Buffer, pH 7.4
Preservatives	0.05% Sodium Azide