

## Product Datasheet

### Phencyclidine (PCP)-BSA Antigen LA018

<b>Product Name</b>	Phencyclidine (PCP)-BSA Antigen
<b>Catalog Number</b>	LA018
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Buffer</b>	Distilled Water
<b>Matched Pairs</b>	HM927 or HM929 (detection)
<b>Preservatives</b>	0.01% Sodium Azide