

Product Datasheet

Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant ? 90% LA026

Product Name	Hepatitis C Virus (HCV) Chimeric Antigen, Recombinant \geq 90%
Catalog Number	LA026
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	\geq 90%
Source	
Buffer	20 mM phosphate buffer, 5 mM EDTA, 0.05% SDS, pH 7.4
Specificity	NS3, Core, NS4
Expression System	E. coli
Preservatives	None