

Product Datasheet

Treponema pallidum (TpN47) Antigen, Recombinant ≥85% LA029

Product Name	Treponema pallidum (TpN47) Antigen, Recombinant ≥85%
Catalog Number	LA029
Tested Applications	ELISA, WB
Alternative Names	Syphilis
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	≥85%
Source	
Buffer	25 mM Tris-HCl, pH 8.0 with 0.1% SDS, 5 mM EDTA
Specificity	TpN47
Expression System	E. coli