

Product Datasheet

Treponema pallidum (TpN47) Antigen, Recombinant ≥85% LA029

Product Name Treponema pallidum (TpN47) Antigen, Recombinant ≥85%

Catalog Number LA029

Tested Applications ELISA, WB

Alternative Names Syphilis

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen

Storage -20°C. Avoid repeated freezing and thawing.

Purity ≥85%

Source

Buffer 25 mM Tris-HCl, pH 8.0 with 0.1% SDS, 5 mM EDTA

Specificity TpN47

Expression System E. coli