

## Product Datasheet

### Hepatitis C Virus (HCV-core) Antigen, Recombinant >98% LA034

<b>Product Name</b>	Hepatitis C Virus (HCV-core) Antigen, Recombinant >98%
<b>Catalog Number</b>	LA034
<b>Tested Applications</b>	WB, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	$\geq 95\%$
<b>Source</b>	
<b>Buffer</b>	25 mM Tris-HCl, 0.1% SDS, 5 mM EDTA, pH 8.0 $\pm$ 0.2
<b>Specificity</b>	Core
<b>Expression System</b>	E. coli
<b>Preservatives</b>	None