

## Product Datasheet

### Human Immunodeficiency Virus 1 (HIV 1) Antigen, Recombinant LA038

<b>Product Name</b>	Human Immunodeficiency Virus 1 (HIV 1) Antigen, Recombinant
<b>Catalog Number</b>	LA038
<b>Tested Applications</b>	LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Source</b>	
<b>Buffer</b>	20 mM Phosphate Buffer, 0.1% SDS, pH 7.4
<b>Specificity</b>	gp120 and gp41
<b>Molecular Weight</b>	43 kDa
<b>Matched Pairs</b>	LA249, LA250 or LA779 (capture)
<b>Expression System</b>	E. coli