

## Product Datasheet

### Treponema pallidum (T. pallidum) Chimeric Antigen, Recombinant >90% LA041

<b>Product Name</b>	Treponema pallidum (T. pallidum) Chimeric Antigen, Recombinant >90%
<b>Catalog Number</b>	LA041
<b>Tested Applications</b>	LF
<b>Alternative Names</b>	Syphilis
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	-20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	≥90%
<b>Source</b>	
<b>Buffer</b>	25 mM Tris-HCl, pH 8.3 ± 0.2 with 0.1% SDS
<b>Specificity</b>	TpN15, TpN17, TpN47
<b>Matched Pairs</b>	LA042 (capture)
<b>Expression System</b>	E. coli