

Product Datasheet

Treponema pallidum (T. pallidum) Chimeric Antigen, Recombinant >80% LA047-1

Product Name	Treponema pallidum (T. pallidum) Chimeric Antigen, Recombinant >80%
Catalog Number	LA042
Tested Applications	LF
Alternative Names	Syphilis
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	-20°C
Purity	>80%
Source	
Buffer	20 mM Phosphate Buffer, pH 7.4 ± 0.2 with 0.15M NaCl
Specificity	TpN15, TpN17, TpN47
Matched Pairs	LA041 (detection)
Expression System	E. coli