

Product Datasheet

Human CD138 ELISA Set 851.620.020

Product Name	Human CD138 ELISA Set
Catalog Number	851.620.020
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	8-256 ng/ml
Sensitivity	4.94 ng/ml
Tested Applications	ELISA
Alternative Names	Syndecan-1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes natural human soluble syndecan-1
Sample Types	Serum, Plasma, Cell culture supernatant