

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-p24) Antigen, Recombinant >95% LA050

Product Name	Human Immunodeficiency Virus 1 (HIV 1-p24) Antigen, Recombinant >95%
Catalog Number	LA050
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	>95%
Buffer	Tris-Saline, pH 8.0
Specificity	p24; Gag
Expression System	P. pastoris