

Product Datasheet

Pepsinogen I (PGI) Antigen, Native >90%
LA072

| | |
|---------------------|---|
| Product Name | Pepsinogen I (PGI) Antigen, Native >90% |
| Catalog Number | LA072 |
| Tested Applications | ELISA, Control, WB |
| Brand | EastCoast Bio |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | Tested negative for HIV, HBV, HCV by currently approved FDA methods. However, material should be treated as potentially infectious. |
| Purity | ≥90% |
| Source | Human Stomach |