

# Product Datasheet

## Human IL-13 ELISA Set 851.630.005

Product Name	Human IL-13 ELISA Set
Catalog Number	851.630.005
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	3.12-100 pg/ml
Sensitivity	1.5 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-13
Sample Types	Serum, Plasma, Cell culture supernatant