

Product Datasheet

Troponin T, Cardiac Antigen, Native $\geq 95\%$ LA149

Product Name	Troponin T, Cardiac Antigen, Native $\geq 95\%$
Catalog Number	LA149
Brand	EastCoast Bio
Form/Appearance	Liquid, Native Antigen
Storage	Short Term (< 1 week): -20°C . Long Term: $\leq -70^{\circ}\text{C}$. Avoid repeated freezing and thawing.
Note	Biohazard Testing - Human tissue donors have tested negative for HBsAg, HIV 1 & 2 antibodies, HIV 1 antigen, HCV antibody, HCV antigen, and RPR Syphilis test. Although material has tested negative, no test can guarantee a product to be non-infectious. All human tissue or fluid should be treated as potentially infectious and good laboratory practices used.
Purity	$\geq 95\%$
Source	Human Cardiac Tissue
Buffer	6M Urea, 250mM NaCl, 50mM Tris-HCl, 15mM β -mercaptoethanol, 1mM EDTA, pH 7.2