

## Product Datasheet

### Troponin T, Cardiac Antigen, Native ≥95% LA149

<b>Product Name</b>	Troponin T, Cardiac Antigen, Native ≥95%
<b>Catalog Number</b>	LA149
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid, Native Antigen
<b>Storage</b>	Short Term (< 1 week): -20°C. Long Term: ≤-70°C. Avoid repeated freezing and thawing.
<b>Note</b>	Biohazard Testing - Human tissue donors have tested negative for HBsAg, HIV 1 & 2 antibodies, HIV 1 antigen, HCV antibody, HCV antigen, and RPR Syphilis test. Although material has tested negative, no test can guarantee a product to be non-infectious. All human tissue or fluid should be treated as potentially infectious and good laboratory practices used.
<b>Purity</b>	≥95%
<b>Source</b>	Human Cardiac Tissue
<b>Buffer</b>	6M Urea, 250mM NaCl, 50mM Tris-HCl, 15mM β-mercaptoethanol, 1mM EDTA, pH 7.2