

# Product Datasheet

## Dengue Virus (DENV-Serotype 1) Antigen, Native LA193

Product Name	Dengue Virus (DENV-Serotype 1) Antigen, Native
Catalog Number	LA193
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	-70°C. Avoid repeated freezing and thawing.
Note	Infectious Dengue Virus was not detectable, but as with all biological reagents, this material should be treated as potentially infectious.
Source	E6 cells/Strain Th-Sman
Buffer	10 mM Phosphate Buffer, pH 7.4 with 140 mM NaCl and containing less than 0.1% CHAPS
Preservatives	None