

# Product Datasheet

## Human CD25 ELISA Set 851.650.005

Product Name	Human CD25 ELISA Set
Catalog Number	851.650.005
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	68.75-2220 pg/ml
Sensitivity	8 pg/ml
Tested Applications	ELISA
Alternative Names	IL-2R
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble receptors
Specificity	Recognizes both natural and recombinant human soluble IL-2R
Sample Types	Serum, Plasma, Cell culture supernatant