

## Product Datasheet

### Dengue Virus (DENV-Serotype 3) Antigen, Native LA195

<b>Product Name</b>	Dengue Virus (DENV-Serotype 3) Antigen, Native
<b>Catalog Number</b>	LA195
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-70°C. Avoid repeated freezing and thawing.
<b>Note</b>	Infectious Dengue Virus was not detectable, but as with all biological reagents, this material should be treated as potentially infectious.
<b>Source</b>	E6 cells/Strain H-87
<b>Buffer</b>	10 mM Phosphate Buffer, pH 7.4 with 140 mM NaCl and containing less than 0.1% CHAPS
<b>Preservatives</b>	None