
Product Datasheet

Human CD25 ELISA Set 851.650.010

Product Name	Human CD25 ELISA Set
Catalog Number	851.650.010
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	68.75-2220 pg/ml
Sensitivity	8 pg/ml
Tested Applications	ELISA
Alternative Names	IL-2R
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble receptors
Specificity	Recognizes both natural and recombinant human soluble IL-2R
Sample Types	Serum, Plasma, Cell culture supernatant