

# Product Datasheet

## Reovirus Serotype 3 Antigen, Native Murine LA202

Product Name	Reovirus Serotype 3 Antigen, Native Murine
Catalog Number	LA202
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	≤-70°C
Note	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
Source	LLC Cells/Strain Abney
Buffer	50 mM Glycine, pH 9 with <0.2% CHAPS and <0.01% SDS
Preservatives	None