

Product Datasheet

Cytomegalovirus (CMV) Antigen, Native  
LA210

Product Name	Cytomegalovirus (CMV) Antigen, Native
Catalog Number	LA210
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	≤-70°C
Note	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
Source	MRC-5 Cells/Strain AD169
Buffer	50 mM Glycine, pH 9 with 140 mM NaCl and <0.1% CHAPS
Preservatives	None