

## Product Datasheet

### C-Reactive Protein, Canine (cCRP) Antigen, Native $\geq 98\%$ LA216

<b>Product Name</b>	C-Reactive Protein, Canine (cCRP) Antigen, Native $\geq 98\%$
<b>Catalog Number</b>	LA216
<b>Tested Applications</b>	WB, ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid, affinity purified
<b>Storage</b>	$\leq -20^{\circ}\text{C}$ . Avoid repeated freezing and thawing.
<b>Purity</b>	$\geq 98\%$ (20 kDa bands with SDS-PAGE)
<b>Source</b>	Dog serum
<b>Buffer</b>	100 mM Tris HCl, pH 7.5 with 200 mM NaCl, and 2 mM $\text{CaCl}_2$
<b>Specificity</b>	Strongly reactive with tested Canine CRP antibodies
<b>Preservatives</b>	0.05% Sodium Azide