

## Product Datasheet

### Adenovirus Type 2, Canine Native Antigen LA224

<b>Product Name</b>	Adenovirus Type 2, Canine Native Antigen
<b>Catalog Number</b>	LA224
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Native Antigen
<b>Storage</b>	≤ -70°C
<b>Note</b>	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
<b>Source</b>	A-72 Cells/Virus Serotype 2
<b>Buffer</b>	50 mM Glycine, pH 9 with 140 mM NaCl and <0.1% CHAPS
<b>Preservatives</b>	None