

Product Datasheet

Adenovirus Type 1, Canine Native Antigen LA225

Product Name	Adenovirus Type 1, Canine Native Antigen
Catalog Number	LA225
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Native Antigen
Storage	$\leq -70^{\circ}\text{C}$
Note	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
Source	A-72 Cells/Virus Serotype 1
Buffer	50 mM Glycine, pH 9 with 140 mM NaCl and <0.1% CHAPS
Preservatives	None