

Product Datasheet

Adenovirus, Mouse Native Antigen LA230

Product Name	Adenovirus, Mouse Native Antigen
Catalog Number	LA230
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Native Antigen
Storage	≤ -70°C
Note	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
Source	NIH 3T3 Cells/Virus Strain FL
Buffer	50 mM Glycine, pH 9 with 140 mM NaCl, <0.3% CHAPS and <0.01% SDS
Preservatives	None