

Product Datasheet

Pneumonia Virus of Mice (PVM) Antigen, Native LA231

Product Name	Pneumonia Virus of Mice (PVM) Antigen, Native
Catalog Number	LA231
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	≤-70°C
Note	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
Source	BHK-21 cells/Strain No. 15
Buffer	50 mM Glycine, pH 9 with <0.2% CHAPS and 0.01% SDS
Preservatives	None