

Product Datasheet

Pneumonia Virus of Mice (PVM) Antigen, Native LA231

| | |
|----------------------------|--|
| Product Name | Pneumonia Virus of Mice (PVM) Antigen, Native |
| Catalog Number | LA231 |
| Tested Applications | ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | Liquid |
| Storage | ≤-70°C |
| Note | Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious. |
| Source | BHK-21 cells/Strain No. 15 |
| Buffer | 50 mM Glycine, pH 9 with <0.2% CHAPS and 0.01% SDS |
| Preservatives | None |