

Product Datasheet

Coronavirus, Canine Antigen, Native LA235

Product Name	Coronavirus, Canine Antigen, Native
Catalog Number	LA235
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	≤ -70°C
Note	Although material has been inactivated, as with all biological reagents, it should be treated as potentially infectious.
Source	A-72 cells/Virus strain 1-71
Buffer	0.05M Glycine, pH 9 with <0.1% CHAPS and <0.01% SDS
Preservatives	None