

Product Datasheet

Human CD106 ELISA Set 851.660.020

Product Name	Human CD106 ELISA Set
Catalog Number	851.660.020
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	1.56-50 ng/ml
Sensitivity	0.6 ng/ml
Tested Applications	ELISA
Alternative Names	VCAM-1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human adhesion molecules
Specificity	Recognizes both natural and recombinant human soluble VCAM-1
Sample Types	Serum, Plasma, Cell culture supernatant