

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1) Chimeric Antigen, Recombinant $\geq 80\%$ LA249

Product Name	Human Immunodeficiency Virus 1 (HIV 1) Chimeric Antigen, Recombinant $\geq 80\%$
Catalog Number	LA249
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	2-8°C
Purity	$\geq 80\%$
Source	
Buffer	Tris-HCl buffer, pH 8.5 \pm 0.2 with 3 mM EDTA
Specificity	gp120 and gp41
Matched Pairs	LA038 (detection)
Expression System	E. coli