

# Product Datasheet

## Human IL-5 ELISA Set 851.670.001

|                     |  |
|---------------------|--|
| Product Name        | Human IL-5 ELISA Set   |
| Catalog Number      | 851.670.001  |
| Description         | Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition. |
| Range               | 7.8-250 pg/ml  |
| Sensitivity         | 5 pg/ml  |
| Tested Applications | ELISA  |
| Alternative Names   | NULL   |
| Brand               | Diaclone   |
| Storage             | Store at +4°C  |
| Cross Reactivity    | No cross reactivity with other human cytokines   |
| Specificity         | Recognizes both natural and recombinant human IL-5   |
| Sample Types        | Serum, Plasma, Cell culture supernatant  |