

# Product Datasheet

## Feline Leukemia Virus (p27) Antigen, Recombinant >90% LA254L

Product Name	Feline Leukemia Virus (p27) Antigen, Recombinant >90%
Catalog Number	LA254L
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Lyophilized Recombinant Antigen (lyophilized from PBS with non-protein stabilizers)
Storage	Lyophilized: 2-8°C. Reconstituted: Short Term 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Reconstitution: Rehydrate with 0.2 ml of reagent grade water
Purity	>90% by SDS-PAGE
Source	
Specificity	amino acids 272-519 with 6xHIS tag
Expression System	E. coli
Preservatives	None