

## Product Datasheet

### Feline Leukemia Virus (p27) Antigen, Recombinant >90% LA254L

<b>Product Name</b>	Feline Leukemia Virus (p27) Antigen, Recombinant >90%
<b>Catalog Number</b>	LA254L
<b>Tested Applications</b>	WB, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Lyophilized Recombinant Antigen (lyophilized from PBS with non-protein stabilizers)
<b>Storage</b>	Lyophilized: 2-8°C. Reconstituted: Short Term 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Reconstitution: Rehydrate with 0.2 ml of reagent grade water
<b>Purity</b>	>90% by SDS-PAGE
<b>Source</b>	
<b>Specificity</b>	amino acids 272-519 with 6xHIS tag
<b>Expression System</b>	E. coli
<b>Preservatives</b>	None