

## Product Datasheet

### Zopiclone-BSA Antigen LA269

<b>Product Name</b>	Zopiclone-BSA Antigen
<b>Catalog Number</b>	LA269
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Buffer</b>	Distilled Water
<b>Matched Pairs</b>	HM310 (detection)
<b>Preservatives</b>	0.01% Sodium Azide