

Product Datasheet

Tricyclic Antidepressant (TCA)-BSA Antigen LA272

Product Name	Tricyclic Antidepressant (TCA)-BSA Antigen
Catalog Number	LA272
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM336 (detection)
Preservatives	0.01% Sodium Azide