

Product Datasheet

Malaria (pf-HRP-2) Antigen, Recombinant $\geq 95\%$ LA278

Product Name	Malaria (pf-HRP-2) Antigen, Recombinant $\geq 95\%$
Catalog Number	LA278
Tested Applications	ELISA, LF, Control
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (His-tag)
Storage	$\leq -20^{\circ}\text{C}$. Avoid repeated freezing and thawing.
Purity	$\geq 95\%$ by SDS-PAGE
Buffer	50 mM Tris, pH 8.0 with 100 mM Sodium Chloride and 0.2M Sucrose
Molecular Weight	~ 75 kDa
Matched Pairs	HM505, HM506 or HM507
Expression System	E. coli