

Product Datasheet

Ketamine-BSA Antigen LA280

Product Name	Ketamine-BSA Antigen
Catalog Number	LA280
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	20 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM390 or HM391 (detection)
Preservatives	0.01% Sodium Azide