

Product Datasheet

Human IL-5 ELISA Set 851.670.020

Product Name	Human IL-5 ELISA Set
Catalog Number	851.670.020
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	7.8-250 pg/ml
Sensitivity	5 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-5
Sample Types	Serum, Plasma, Cell culture supernatant