
Product Datasheet

EDDP-BSA Antigen LA282

Product Name	EDDP-BSA Antigen
Catalog Number	LA282
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM381 (detection)
Preservatives	0.01% Sodium Azide