

Product Datasheet

Buprenorphine (BUP)-BSA Antigen LA283

Product Name	Buprenorphine (BUP)-BSA Antigen
Catalog Number	LA283
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	20 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM376 or HM377 (detection)
Preservatives	0.01% Sodium Azide