

Product Datasheet

Dehydroepiandrosterone sulfate (DHEA-S)-BSA Antigen LA285

| | |
|----------------------------|---|
| Product Name | Dehydroepiandrosterone sulfate (DHEA-S)-BSA Antigen |
| Catalog Number | LA285 |
| Tested Applications | ELISA, LF |
| Brand | EastCoast Bio |
| Form/Appearance | Liquid; Antigen-BSA conjugate |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Matched Pairs | HM378, HM379 or HM380 (detection) |
| Preservatives | 0.05% Proclin 950 |