

Product Datasheet

Pregnanediol Glucuronide (PDG)-BSA Antigen LA286

Product Name	Pregnanediol Glucuronide (PDG)-BSA Antigen
Catalog Number	LA286
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	≥ 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, 0.05% sodium azide, pH 7.2 ± 0.2
Matched Pairs	HM386 or HM387 (detection)