

Product Datasheet

Malaria (pf-HRP-2) Antigen, Recombinant >90% LA287

Product Name	Malaria (pf-HRP-2) Antigen, Recombinant >90%
Catalog Number	LA287
Tested Applications	ELISA, LF, Control
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term (\leq 1 week): 2-8°C. Long Term: -70°C. Avoid repeated freezing and thawing.
Purity	\geq 90%
Buffer	20 mM Tris, pH 8.0 with 20% glycerol
Molecular Weight	46 kDa
Expression System	E. coli
Preservatives	0.05% Proclin 300