

Product Datasheet

Human IL-7 ELISA Set 851.680

Product Name	Human IL-7 ELISA Set
Catalog Number	851.680
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	6.25-200 pg/ml
Sensitivity	3 pg/ml
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reaction with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-7
Sample Types	Serum, Plasma, Cell culture supernatant
References	Kraus, H. et al., J. Immunol.,2014;192(3):1044-1054 - Pubmed link Francois,B. et al.,JCI Insight. 2018 Mar 8;3(5). pii: 98960. - Pubmed link